

LSRC85-5

4.6 Cubic Ft. ULT Chest Freezer



-86°C Temperature Chest Freezer

Temperature Range -40°C to -86°C

Exterior w x d x h
28.5" x 25.8" x 34.1"


Internal w x d x h
20.5x17.7x25.6

Capacity
4.6 cu. t.

Weight (lb)
124 lbs

Used by universities, hospitals and research institutes for long term storing of blood, skin cells, DNA, organs, and vaccines

Digital controller with potential free contact (which makes it possible to connect it to an external alarm)



- Gross volume, cu.ft 4.6
- Refrigerant: Mix of Hydro Carbons
- Energy consumption / 24 h kwh at 10 v / 60 hz and 40 % running time: 4.6
- Input power (watt) at 110 volt / 60 hz: 480
- External dimensions (WxDxH) 28.5x25.8x34.1
- Internal dimensions (WxDxH) 20.5x17.7x25.6
- Type of defrost: Manual
- Power: 110V/60Hz
- Net weight, lb: 124
- Gross weight, lb: 139
- Insulation, mm: 3.9
- Outer cabinet: White



This unit comes with a standard 2 year warranty.

LSR cold storage solutions

(800) 506-7470 | lsrfreezer.com | info@lsrfreezer.com