



► Main Features

- 10X Wide Field Monocular Dual Head.
- LED lighting with intensity adjustment.
- Compact, lightweight and easy to transport.
- Stage with clamps, stop with height adjustment, vernier and millimeter scale.
- Colored objective rings for easy identification.
- Non-slip lenses with metallic finish.
- Finished in epoxy paint.
- Focus with tension adjustment.
- Head inclined at 45° and rotatable 360°.

► Applications

- Education

► Accessories Included

- AC/CD Adapter
- Dust cover
- Immersion oil
- Color filters green, yellow and blue

► TECHNICAL DATA

Model	VE-M5DTH
Eyepiece	WF 10X/18 mm with block screw and pointer.
Nosepiece	Quadruple
Objectives	4X, 10X, 40X y 100X (oil) (S)
Optical System	Achromatic
Stage	Double layer with coaxial movement X-Y 130 x 130 mm / 5.1" x 5.1".
Condenser	Abbe, N.A 1.25
Diaphragm	Iris type with filter holder.
Focus	Coarse and fine anti-slip control knob.
Illumination	LED.
Dimensions	145 x 350 x 230 mm/ 5.7" x 13.7 x 9.0".
Power Supply	Three 1.2V AA rechargeable batteries (included). 110/240, 50/60 Hz AC Adapter.
Warranty	10 years