

# LSRP-RG-72

60 Cubic Ft. Glass Door Refrigerator



## Triple Glass Door Refrigerator

### Accurate Control System

This 2°C~8°C medical/vaccine fridge comes equipped with a high precision temperature control system and finely calibrated sensors to ensure temperature uniformity within  $\pm 2^{\circ}\text{C}$ . We designed this refrigerator with a high brightness digital temperature and humidity display for easy to read automatic temperature control and display accuracy within 0.1°C.

### Powerful Refrigeration System

This model is equipped with a multi-duct vortex refrigeration system with a finned evaporator which can prevent frosting and improves temperature uniformity to a large extent. The air cooling system with auto-defrost utilizes HCFC-FREE refrigerant to deliver more effective refrigeration and ensures environmentally friendly operation.

### Human-Oriented Design

The inside of the refrigerator is designed with a bright LED lighting system for easy visibility within the unit. The light will turn on while door is opened, and the light will go off once the door is closed. The exterior cabinet is made with high quality powder coated steel, and the interior stainless steel is easy to clean and maintain.

### Temperature Range(°C)

2° to 8°

### Exterior w x d x h (in)

70.7 x 34.6 x 78.3

### Internal w x d x h (in)

66.3 x 26.4 x 59.6

### Capacity (Cu.Ft)

60

### Weight (lbs)

710

## Model

Cabinet Type	Upright
Capacity (Cu.Ft)	60
Internal Size (WxDxH) (inches)	66.3x26.4x59.6
External Size (WxDxH) (inches)	70.7x34.6x78.3
Package Size (WxDxH) (inches)	76.1x38.9x85.4
NW/GW(Lbs)	710/950

## Performance

Temperature Range (C°)	2°~ 8°
Ambient Temperature (C°)	16° - 32°
Controller	Microprocessor
Display	Digital display

## Refrigeration

Compressor	Single
Cooling Method	Air cooling
Defrost Mode	Automatic
Refrigerant	R507
Insulation Thickness(mm)	55

## Construction

External Material	Powder Coated Steel
Inner Material	304 Stainless Steel
Shelves	18 (Coated steel wire)
Door Lock with Key	Yes
Lighting	LED
Access Port	1pc. Ø 25 mm
Casters	6 (3 casters with brake)
Data Logging	USB/Record every 10 minute / 2 years

## Alarm

Temperature	High/Low temperature, High ambient temperature,
Electrical	Power failure, and Low battery
System	Sensor error, Door ajar, Built-in datalogger USB Failure, Main board communication error, Remote alarm

## Electrical

Power Supply(V/HZ)	110/60
Rated Current(A)	13.21

Due to periodic changes in designs, methods, procedures, policies and regulations, the specifications contained in this sheet are subject to change without notice. While we exercise good faith efforts to provide information that is accurate, we are not responsible for errors or omissions in information provided or conclusions reached as a result of using the specifications. By using the information provided, the user assumes all risks in connection with such use.

